

LCP GREENS MEDIAL PROXIMAL
TIBIA -C PLATE 3.5MM



LCP GREENS MEDIAL PROXIMAL TIBIA -C PLATE 3.5MM		
CAT. NO.	PRODUCT DESCRIPTION	MAT.
632.030L	LCP GREENS MEDIAL PROXIMAL TIBIA -C PLATE 3.5MM LEFT	TT
632.030R	LCP GREENS MEDIAL PROXIMAL TIBIA -C PLATE 3.5MM RIGHT	TT

Specification

- Product is available in TT.

Indication for use

• The LCP Greens Medial Proximal Tibia-C Plate 3.5 mm is indicated for fixation of fractures, osteotomies, and nonunions of the medial proximal tibia. It is specifically designed for medial tibial plateau and metaphyseal fracture stabilization. The locking compression design provides angular stability and supports fixation in both normal and osteoporotic bone. It is suitable for trauma and reconstructive orthopaedic procedures.

