



CAT. NO.	METAL	STEP SPINALROD 5.5 /3.5 mm
665.810	TT	Step Spinal Rod - Ø 5.5 /3.5 mm x 100 mm
665.812	TT	Step Spinal Rod - Ø 5.5 /3.5 mm x 120 mm
665.814	TT	Step Spinal Rod - Ø 5.5 /3.5 mm x 140 mm
665.815	TT	Step Spinal Rod - Ø 5.5 /3.5 mm x 150 mm
665.816	TT	Step Spinal Rod - Ø 5.5 /3.5 mm x 160 mm
665.818	TT	Step Spinal Rod - Ø 5.5 /3.5 mm x 180 mm
665.820	TT	Step Spinal Rod - Ø 5.5 /3.5 mm x 200 mm
665.825	TT	Step Spinal Rod - Ø 5.5 /3.5 mm x 250 mm
665.827	TT	Step Spinal Rod - Ø 5.5 /3.5 mm x 270 mm
665.830	TT	Step Spinal Rod - Ø 5.5 /3.5 mm x 300 mm
665.834	TT	Step Spinal Rod - Ø 5.5 /3.5 mm x 340 mm
665.835	TT	Step Spinal Rod - Ø 5.5 /3.5 mm x 350 mm
665.840	TT	Step Spinal Rod - Ø 5.5 /3.5 mm x 400 mm
665.845	TT	Step Spinal Rod - Ø 5.5 /3.5 mm x 450 mm
665.850	TT	Step Spinal Rod - Ø 5.5 /3.5 mm x 500 mm

Specification

- Spinal Rod available in Titanium.
- Available length 100 & 500mm with 5.5/3.5mm diameter.

Indication for use

- It is a component of Spine Fixation System intended for use in the pedicle Screw fixation or in the anterolateral lateral spine fixation system.

