



ANGLE PLATE OF OSTEOTOMY 3.5 MM

SR. NO.	PRODUCT DESCRIPTION	MAT.
1	ANGLE PLATE OF OSTEOTOMY 3.5MM 4 HOLES X 20 MM 90°	TT
2	ANGLE PLATE OF OSTEOTOMY 3.5MM 4 HOLES X 25 MM 90°	TT
3	ANGLE PLATE OF OSTEOTOMY 3.5MM 4 HOLES X 30 MM 90°	TT
4	ANGLE PLATE OF OSTEOTOMY 3.5MM 3 HOLES X 30 MM 120°	TT
5	ANGLE PLATE OF OSTEOTOMY 3.5MM 3 HOLES X 35 MM 120°	TT

Specification

- Product is available in Titanium material.

Indication for use

• The Angle Plate for Osteotomy (3.5 MM) is indicated for fixation and stabilization of bone segments following corrective osteotomy procedures. It is designed to maintain proper angular alignment and support bone healing during reconstruction or deformity correction. The plate provides rigid internal fixation for various orthopaedic applications involving long bones. It is suitable for trauma, corrective, and reconstructive surgical procedures.