



## DHS PLATE 3.5MM

SR. NO.	PRODUCT DESCRIPTION	MAT.
1	DHS PLATE 3.5 FOR 3 -7 YEARS INFANT - 3 HOLES X 48 MM 115°	TT
2	DHS PLATE 3.5 FOR 3 -7 YEARS INFANT - 3 HOLES X 48 MM 125°	TT
3	DHS PLATE 3.5 FOR 3 -7 YEARS INFANT - 3 HOLES X 48 MM 135°	TT
4	DHS PLATE 3.5 FOR 3 -7 YEARS INFANT - 3 HOLES X 48 MM 145°	TT

### Specification

- Product is available in Titanium material.

### Indication for use

- The DHS Plate 3.5 mm for Infant is indicated for fixation of proximal femoral fractures and osteotomies in infant and young paediatric patients. It is designed to provide stable internal fixation while accommodating the smaller anatomy of growing bones. The plate supports proper alignment, stabilization, and healing of the proximal femur. It is suitable for paediatric trauma and corrective orthopaedic procedures.