



GROWTH PLATE		
SR. NO.	PRODUCT DESCRIPTION	MAT.
1	GROWTH PLATE 12 MM X 3.5	TT
2	GROWTH PLATE 16 MM X 3.5	TT
3	GROWTH PLATE 16 MM X 3.5	TT
4	GROWTH PLATE 16 MM X 3.5	TT
5	GROWTH PLATE 16 MM X 3.5 (4 HOLES)	TT
6	GROWTH PLATE 12 MM X 4.5	TT
7	GROWTH PLATE 16 MM X 4.5	TT
8	GROWTH PLATE 16 MM X 4.5	TT
9	GROWTH PLATE 16 MM X 4.5	TT
10	GROWTH PLATE 16 MM X 4.5 (4 HOLES)	TT

**Specification**

- Product is available in Titanium material.

**Indication for use**

- The Growth Plate is indicated for guided growth procedures in paediatric orthopaedic patients to correct angular deformities of long bones. It provides temporary stabilization across the growth plate to allow controlled growth modulation and gradual deformity correction. The implant is designed to maintain proper alignment while preserving natural bone growth. It is suitable for paediatric lower limb deformity correction procedures.

