

LCP DISTAL FEMUR PLATE FOR SUPRACONDYLAR OSTERFOMY



| LCP DISTAL FEMUR PLATE FOR SUPRACONDYLAR OSTERFOMY | | |
|--|--|------|
| SR. NO. | PRODUCT DESCRIPTION | MAT. |
| 1 | LCP DISTAL FEMUR PLATE FOR SUPRACONDYLAR OSTERFOMY 3 HOLES 90° X 3.5 | TT |
| 2 | LCP DISTAL FEMUR PLATE FOR SUPRACONDYLAR OSTERFOMY 5 HOLES 90° X 3.5 | TT |
| 3 | LCP DISTAL FEMUR PLATE FOR SUPRACONDYLAR OSTERFOMY 7 HOLES 90° X 3.5 | TT |
| 4 | LCP DISTAL FEMUR PLATE FOR SUPRACONDYLAR OSTERFOMY 3 HOLES 90° X 5.0 | TT |
| 5 | LCP DISTAL FEMUR PLATE FOR SUPRACONDYLAR OSTERFOMY 5 HOLES 90° X 5.0 | TT |
| 6 | LCP DISTAL FEMUR PLATE FOR SUPRACONDYLAR OSTERFOMY 7 HOLES 90° X 5.0 | TT |

Specification

- Product is available in Titanium material.

Indication for use

• The LCP Distal Femur Plate for Supracondylar Osterfomy is indicated for fixation of distal femoral osteotomies and supracondylar femur fractures. It is designed to provide stable internal fixation and angular stability during corrective alignment procedures. The locking compression technology supports secure fixation in both normal and osteoporotic bone. It is suitable for deformity correction, trauma management, and reconstructive orthopaedic surgeries.