

# LCP PROXIMAL FEMORAL PLATE FOR VARUS OSTEOFOMY



LCP PROXIMAL FEMORAL PLATE FOR VARUS OSTEOFOMY		
SR. NO.	PRODUCT DESCRIPTION	MAT.
1	LCP PROXIMAL FEMORAL PLATE FOR VARUS OSTEOFOMY 3 HOLES 100° X 2.7	TT
2	LCP PROXIMAL FEMORAL PLATE FOR VARUS OSTEOFOMY 3 HOLES 110° X 2.7	TT
3	LCP PROXIMAL FEMORAL PLATE FOR VARUS OSTEOFOMY 3 HOLES 100° X 3.5	TT
4	LCP PROXIMAL FEMORAL PLATE FOR VARUS OSTEOFOMY 3 HOLES 110° X 3.5	TT
5	LCP PROXIMAL FEMORAL PLATE FOR VARUS OSTEOFOMY 3 HOLES 100° X 5.0	TT
6	LCP PROXIMAL FEMORAL PLATE FOR VARUS OSTEOFOMY 3 HOLES 110° X 5.0	TT

## Specification

- Product is available in Titanium material.

## Indication for use

- The LCP Proximal Femoral Plate for Varus Osteotomy – 3 Holes is indicated for fixation of proximal femoral osteotomies in paediatric patients. It is specifically designed for corrective varus osteotomy procedures to treat hip deformities and paediatric orthopaedic conditions. The locking compression design provides stable fixation and angular stability in growing bone. It supports proper alignment, stabilization, and healing of the proximal femur.