


2.7MM MANDIBLE EMERGENCY CORTEX SCREW SELF TAPPING



2.7MM MANDIBLE EMERGENCY CORTEX SCREW SELF TAPPING	
CAT NO.	PRODUCT DESCRIPTION
706.405	2.7MM MANDIBLE EMERGENCY CORTEX SCREW SELF TAPPING 5.0 MM
706.406	2.7MM MANDIBLE EMERGENCY CORTEX SCREW SELF TAPPING 6.0 MM
706.408	2.7MM MANDIBLE EMERGENCY CORTEX SCREW SELF TAPPING 8.0 MM
706.410	2.7MM MANDIBLE EMERGENCY CORTEX SCREW SELF TAPPING 10 MM
706.412	2.7MM MANDIBLE EMERGENCY CORTEX SCREW SELF TAPPING 12 MM
706.414	2.7MM MANDIBLE EMERGENCY CORTEX SCREW SELF TAPPING 14 MM
706.416	2.7MM MANDIBLE EMERGENCY CORTEX SCREW SELF TAPPING 16 MM
706.418	2.7MM MANDIBLE EMERGENCY CORTEX SCREW SELF TAPPING 18 MM

Specification

- Screw is available in titanium.

Indication for use

- A Mandible Emergency Cortex Screw is a self-tapping cortical screw used in oral and maxillofacial emergencies, primarily for the immediate fixation of mandibular fractures or to secure plates during trauma surgeries.