



CAT. NO.	L-PLATE 2+3 HOLES
706.506L	L-PLATE 2+3 HOLES LEFT THICKNESS 0.6MM
706.506R	L-PLATE 2+3 HOLES RIGHT THICKNESS 0.6MM
707.006L	L-PLATE 2+3 HOLES LEFT THICKNESS 0.7MM
707.006R	L-PLATE 2+3 HOLES RIGHT THICKNESS 0.7MM
707.606L	L-PLATE 2+3 HOLES LEFT THICKNESS 0.8MM
707.606R	L-PLATE 2+3 HOLES RIGHT THICKNESS 0.8MM

## Specification

- Plate is available in Titanium.

## Indication for use

- The L-Plate (2+3 holes) is used in oral and maxillofacial surgery, particularly for fixation of facial bone fractures, where a right-angle configuration is beneficial. It is commonly used for mandibular angle fractures, orbital rim fractures, or zygomaticomaxillary complex (ZMC) repairs.