



CAT NO.	LCP MANDIBLE RECONSTRUCTION PLATE, STRAIGHT
709.812	LCP MANDIBLE Reconstruction Plate, straight, 12 holes THICKNESS 2.5MM
709.820	LCP MANDIBLE Reconstruction Plate, straight, 20 holes THICKNESS 2.5MM
710.512	LCP MANDIBLE Reconstruction Plate 12 holes, Thickness 2.8 MM
710.520	LCP MANDIBLE Reconstruction Plate, straight, 20 holes, Thickness 2.8 MM

Specification

- Plate is available in Titanium.

Indication for use

- The straight LCP Mandible Reconstruction Plate is designed for the internal fixation of mandibular fractures or defects, particularly in the horizontal ramus and interdental spaces. Its straight configuration allows for stable fixation with minimal bone contact, promoting optimal healing.