

MAXILLOFACIAL EMERGENCY SCREW SELF TAPPING 2.3 MM



MAXILLOFACIAL EMERGENCY SCREW SELF TAPPING 2.3 MM

CAT NO.	PRODUCT DESCRIPTION
706.283	MAXILLOFACIAL EMERGENCY SCREW SELF TAPPING 2.3MM X 3 MM
706.284	MAXILLOFACIAL EMERGENCY SCREW SELF TAPPING 2.3MM X 4 MM
706.285	MAXILLOFACIAL EMERGENCY SCREW SELF TAPPING 2.3MM X 5 MM
706.286	MAXILLOFACIAL EMERGENCY SCREW SELF TAPPING 2.3MM X 6 MM
706.287	MAXILLOFACIAL EMERGENCY SCREW SELF TAPPING 2.3MM X 7 MM
706.288	MAXILLOFACIAL EMERGENCY SCREW SELF TAPPING 2.3MM X 8 MM
706.289	MAXILLOFACIAL EMERGENCY SCREW SELF TAPPING 2.3MM X 9 MM
706.290	MAXILLOFACIAL EMERGENCY SCREW SELF TAPPING 2.3MM X 10 MM
706.291	MAXILLOFACIAL EMERGENCY SCREW SELF TAPPING 2.3MM X 11 MM
706.292	MAXILLOFACIAL EMERGENCY SCREW SELF TAPPING 2.3MM X 12 MM
706.294	MAXILLOFACIAL EMERGENCY SCREW SELF TAPPING 2.3MM X 14 MM
706.293	MAXILLOFACIAL EMERGENCY SCREW SELF TAPPING 2.3MM X 16 MM
706.298	MAXILLOFACIAL EMERGENCY SCREW SELF TAPPING 2.3MM X 18 MM

Specification

- Screw is available in titanium.

Indication for use

- The Maxillofacial Emergency Screw Self-Tapping is indicated for temporary or urgent fixation of facial bone fractures, osteotomies, or trauma cases. Its self-tapping design enables quick and secure placement without pre-drilling, ideal for emergency surgical interventions.

