



MESH PLATE FOR 1.5 SCREWS

CAT. NO.	PRODUCT DESCRIPTION
704.901	MESH PLATE FOR 1.5 SCREWS 0.6MM X 25MM X 25MM
704.902	MESH PLATE FOR 1.5 SCREWS 0.6MM X 25MM X 50MM
704.920	MESH PLATE FOR 1.5 SCREWS 0.6MM X 100MM X 100MM

Specification

- Product is available in TT.

Indication for use

- The Mesh Plate for 1.5 Screws is used in maxillofacial and cranial procedures for the fixation and stabilization of small bone fragments. It allows contouring to complex anatomical surfaces and is secured with 1.5 mm screws to maintain alignment and support bone healing.