



## HERBERT SCREW SYSTEM

SR. NO.	PRODUCT DESCRIPTION	MAT.
1	3.5MM CANNULATED HERBERT SCREW	TT
2	2.4MM HEADLESS COMPRESSION SCREW CANNULATED (HERBERT)	TT
3	4.0MM HEADLESS COMPRESSION SCREW CANNULATED(HERBERT)	TT



HERBERT SCREW INSTRUMENT SET		
SR. NO.	NAME OF INSTRUMENTS	QTY
1	CANNULATED SCREW DRIVER 1.5MM	1
2	CANNULATED SCREW DRIVER 2.0MM	1
3	CANNULATED SCREW DRIVER 2.5MM	1
4	DEPTH GAUGE	1
5	GUIDE WIRE, PLAIN- Ø 0.8 MM	2
6	GUIDE WIRE, PLAIN- Ø 1.0 MM	2
7	GUIDE WIRE, PLAIN- Ø 1.2 MM	2
8	CANNULATED DRILL TAP 2.0MM	1
9	CANNULATED DRILL TAP 2.5MM	1
10	CANNULATED DRILL TAP 3.5MM	1
11	CANNULATED DRILL TAP 4.5MM	1
12	CANNULATED DRILL BIT 1.8MM	1
13	CANNULATED DRILL BIT 2.0MM	1
14	CANNULATED DRILL BIT 2.7MM	1
15	CANNULATED DRILL BIT 3.2MM	1
16	DRILL SLEEVE	1
17	PROTECTION SLEEVE	1
18	JIG FOR HERBURT	1
19	HERBERT SCREW BOX	1

## Specification

- All implants are available in TT and SS.

## Indication for use

- The Herbert screw instrument set is indicated for internal fixation of small bone fractures, particularly scaphoid fractures, nonunions, and osteotomies. It provides stable, headless compression fixation ideal for areas requiring minimal soft tissue disruption.

