


HTO PLATE INSTRUMENT SET	
SR. NO.	NAME OF INSTRUMENTS
1	ORTHOPEDIC DISTRACTION FORCEPS
2	HEIGHT MEASURING DEVICE
3	ANGLE DISTRACTION DEVICE
4	OSTEOTOMY CHISEL 10MM
5	OSTEOTOMY CHISEL 15MM
6	OSTEOTOMY CHISEL 20MM
7	OSTEOTOMY CHISEL 25MM
8	GUIDE PIN Φ 2.5
9	GUIDE PIN Φ 2
10	DRILL GUIDE
11	BONE DRILL
12	BONE DRILL
13	DRILL STOPPER Φ 4.1
14	L-SHAPED HEXAGONAL WRENCH
15	DRILL GUIDE
16	TAP Φ 5
17	SCREWDRIVER
18	BONE DRILL Φ 3.2
19	TAP HA4.5
20	QUICK-CHANGE HANDLE
21	DRILL GUIDE
22	HEXAGONAL WRENCH SW3.5
23	4.0N.M TORQUE HANDLE
24	DEPTH GAUGE
25	ADJUSTABLE FIXATION BRACE (A)
26	ADJUSTABLE FIXATION BRACE (B)
27	MEDICAL RULER
28	MEDICAL FIXATION BELT
29	T25 PLUM WRENCH





Specification

- All instruments are available in SS.

Indication for use

- HTO Plate Instrument Set is used in High Tibial Osteotomy procedures to correct knee alignment by offloading the medial compartment. It includes osteotomes, spreaders, drill guides, and locking plates with screws for precise bone cutting and fixation.