



## LARGE CC SCREW IMPLANT SYSTEM

SR. NO.	PRODUCT NAME
1	LARGE CANNULATED SCREWS Ø 6.5MM (F.TH.)
2	LARGE CANNULATED SCREWS Ø 6.5MM/16 MM TH.
3	LARGE CANNULATED SCREWS Ø 6.5MM - 32TH TH.
4	LARGE CANNULATED SCREWS Ø 7.3MM - 16TH TH.
5	LARGE CANNULATED SCREWS Ø 7.3MM - 32TH TH.
6	LARGE CANNULATED SCREWS Ø 7.3MM - FULL TH.

LARGE CC SCREW IMPLANT INSTRUMENT SET		
SR. NO.	NAME OF INSTRUMENT	QTY
1	CLEAN GUIDE PIN $\Phi$ 2.5*250	1
2	SLEEVE FOR GUIDE PIN $\Phi$ 2.5	1
3	SLEEVE FOR DRILL BIT $\Phi$ 5	1
4	DILATOR $\Phi$ 2.5	1
5	TWO-HOLE GUIDER $\Phi$ 2.5	1
6	MULTI-HOLE GUIDER	1
7	CANNULATED DRILL BIT $\Phi$ 5	1
8	CANNULATED TAP $\Phi$ 7.3	1
9	CANNULATED COUNTERSINK DRILL	1
10	CANNULATED SLEEVE SW6	1
11	SKID TAKE-OFF SW4	1
12	LIMITED WRENCH SW3	1
13	CANNULATED SCREWDRIVER HEXSW4.0	1
14	SCREWDRIVER HEX SW4.0	1
15	DEPTH GAUGE FOR GUIDE PIN $\Phi$ 2.5	1
16	GUIDE PIN WITH THREADS 2.5*250	5
17	CASE	1

## Specification

- All implants are available in TT and All instruments are available in SS.

## Indication for use

- The 6.5/7.3 mm Large Cannulated Screw System is intended for fixation of fractures, osteotomies, and arthrodesis of large bones. It is commonly used for femoral neck fractures, tibial plateau fractures, sacroiliac joint fixation, and pelvic fractures. The cannulated design enables accurate placement over a guide wire, ensuring stable fixation with minimal soft tissue damage.

