


THR CEMENTED & UN-CEMENTED ACETABULAR PREPARATION INSTRUMENT SET



THR CEMENTED & UN -CEMENTED ACETABULAR PREPARATION INSTRUMENT SET		
CAT. NO.	NAME OF INSTRUMENTS	QTY
1571.340	ACETABULAR REAMER 40 MM	1
1571.342	ACETABULAR REAMER 42 MM	1
1571.344	ACETABULAR REAMER 44 MM	1
1571.346	ACETABULAR REAMER 46 MM	1
1571.348	ACETABULAR REAMER 48 MM	1
1571.350	ACETABULAR REAMER 50 MM	1
1571.352	ACETABULAR REAMER 52 MM	1
1571.354	ACETABULAR REAMER 54 MM	1
1571.356	ACETABULAR REAMER 56 MM	1
1571.358	ACETABULAR REAMER 58 MM	1
1571.360	ACETABULAR REAMER 60 MM	1
1571.375	TRIAL HEAD 32+5	1
1571.412	TRIAL HEAD 36+12	1
1571.408	TRIAL HEAD 36+8.5	1
1571.405	TRIAL HEAD 36+5	1

THR CEMENTED & UN-CEMENTED ACETABULAR PREPARATION INSTRUMENT SET



1571.395	TRIAL HEAD 36+1.5	1
1571.383	TRIAL HEAD 32+13	1
1571.379	TRIAL HEAD 32+9	1
1571.371	TRIAL HEAD 32+1	1
1571.444	UN CEMENTED TRIAL SHELL 44	1
1571.446	UN CEMENTED TRIAL SHELL 46	1
1571.448	UN CEMENTED TRIAL SHELL 48	1
1571.450	UN CEMENTED TRIAL SHELL 50	1
1571.452	UN CEMENTED TRIAL SHELL 52	1
1571.454	UN CEMENTED TRIAL SHELL 54	1
1571.456	UN CEMENTED TRIAL SHELL 56	1
1571.458	UN CEMENTED TRIAL SHELL 58	1
1571.460	UN CEMENTED TRIAL SHELL 60	1
1571.484	UN CEMENTED TRIAL LINER 44	1
1571.486	UN CEMENTED TRIAL LINER 46	1
1571.488	UN CEMENTED TRIAL LINER 48	1
1571.490	UN CEMENTED TRIAL LINER 50	1
1571.492	UN CEMENTED TRIAL LINER 52	1
1571.494	UN CEMENTED TRIAL LINER 54	1
1571.496	UN CEMENTED TRIAL LINER 56	1
1571.498	UN CEMENTED TRIAL LINER 58	1
1571.500	UN CEMENTED TRIAL LINER 60	1
1571.533	SCREW FOR UN CEMENTED LINER AND CUP	9
1571.534	CLIP WASHER	9
1571.540	CEMENTED ACETABULAR TRAIL CUP 40 MM	1
1571.543	CEMENTED ACETABULAR TRAIL CUP 43 MM	1
1571.547	CEMENTED ACETABULAR TRAIL CUP 47 MM	1
1571.550	CEMENTED ACETABULAR TRAIL CUP 50 MM	1
1571.553	CEMENTED ACETABULAR TRAIL CUP 53 MM	1
1571.436	UNIVERSAL SCREW DRIVER 4.5 MM (FLEXIBLE)	1



THR CEMENTED & UN-CEMENTED ACETABULAR PREPARATION INSTRUMENT SET



900.938	CEMENTED CUP HOLDER	1
1571.432	SOFT DRILL 29 MM	1
1571.430	SOFT DRILL 35 MM	1
1571.438	GREENS DEEP PROBE (DEPTH GAUGE)	1
1571.426	HEXAGONAL SCREW DRIVER 4.5MM	1
1571.424	SCREW HOLDING SPANNER	1
1571.422	ADAPTER FOR REAMER (REAMER HANDLE)	1
1571.210	UN CEMENTED TRAIL CUP INTRODUCER	1
	BOX FOR THR CEMENTED / UN -CEMENTED ACETABULAR PREPARATION	1

Specification

- Instruments are available in SS.

Indication for use

- The THR Cemented & Un-cemented Acetabular Preparation Instrument Set is indicated for use in total hip replacement procedures to prepare the acetabulum for both cemented and un-cemented cup implantation. It facilitates precise reaming, sizing, and contouring of the acetabular cavity. The instruments ensure optimal fit, alignment, and stability of the implant. It is suitable for use in both primary and revision hip arthroplasty surgeries.

