



THR STEM PREPARATION INSTRUMENT SET		
CAT. NO.	NAME OF INSTRUMENTS	QTY
1571.006	STEMP RASP 6.0	1
1571.007	STEMP RASP 7.0	1
1571.008	STEMP RASP 8.0	1
1571.009	STEMP RASP 9.0	1
1571.010	STEMP RASP 10.0	1
1571.011	STEMP RASP 11.0	1
1571.012	STEMP RASP 12.0	1
1571.013	STEMP RASP 13.0	1
1571.014	STEMP RASP 14.0	1
1571.026	TRIAL NECK	1
900.807	QUICK COUPLING TRIAL HANDLE	1

1571.027	CHARNLEY RETRACTOR	1 SET
1571.022	IMPACTOR EXTRACTOR FOR RASP (HANDLE)	1
1571.025	FURCATION HAMMER	1
1571.021	ACETABULAR CURETTE	1
1571.023	STEM IMPACTOR - A	1
1571.024	STEM IMPACTOR	1
	BONE HOOK	1
1571.030	CANAL OPENER (BOX CHISEL)	1
1571.037	LATERALIZER FOR T.H.R .	1
1571.031	CHARNLEY REAMER (MEDULLARY CAVITY RASP)	1
1571.029	MURPHY SKID	1
1571.028	DISTAL PLUNGE LOCATOR (CEMENT RESTRACTOR INTRODUCER)	1
1571.038	OSTEOTOMY GUIDE FOR T.H.R .	1
	BOX FOR T.H.R. STEM PREPARATION INSTRUMENTS	1

Specification

- All instruments are available in SS.

Indication for use

- The THR Stem Preparation Instrument Set is used during Total Hip Replacement (THR) procedures to prepare the femoral canal for implantation of the femoral stem prosthesis. It assists surgeons in reaming, broaching, and shaping the femoral canal to achieve proper alignment, sizing, and stable fixation of the hip stem implant.

