



| WIRE INSTRUMENT SET | | |
|---------------------|--|-----|
| SR. NO. | NAME OF INSTRUMENT | QTY |
| 1 | REDUCTION FORCEPS POINTED - 150MM (SPEED LOCK) | 2 |
| 2 | NOSE PLIER CUM WIRE BENDER CUM WIRE CUTTER - 175MM | 2 |
| 3 | WIRE TIGHTENER | 1 |
| 4 | WIRE CUTTER SMALL | 1 |
| 5 | GUIDE WIRE / K. WIRE INTRODUCER | 1 |
| 6 | JUMBO WIRE CUTTER | 1 |
| 7 | K. WIRE BENDER | 1 |
| 8 | WIRE CUTTER SMALL | 1 |
| 9 | WIRE CUTTER LARGE | 1 |
| 10 | WIRE SLEEVE Ø 2 MM X 2.5 MM | 1 |
| 11 | GUIDE WIRE, THREADED - Ø 2.0 MM | 2 |
| 12 | GUIDE WIRE, THREADED - Ø 3.0 MM | 2 |

Specification

- Product is available in SS.

Indication for use

- The Wire Instrument Set is used in orthopaedic procedures for the insertion, tensioning, and fixation of surgical wires such as K-wires and cerclage wires. It assists in stabilizing fractures and securing bone fragments during fixation procedures.

